

The crowning in level of detail



Still one of the fastest A-silicones

AFFINIS fast, with a maximum total setting time of 2 min 10 sec guarantees optimum efficiency and saves time.



Precise moulding of every detail even under the most difficult conditions

The highly active surfactant system is characterized by excellent surface affinity in a moist environment.



Confidence during impression taking

Self-contouring consistency prevents air bubbles and distortions.















Easy monitoring for success

AFFINIS PRECIOUS wash materials with impressive gold and silver pigmentation for excellent readability and detail reproduction.







AFFINIS® Family

Impression Material

Art. Nr.	Product	Consistency	Processing Time			Total setting time	Handmix	System 50/75	System 360
			Max. processing time	Intraoral setting time					
CW-6501	light body		01:00	02:00		03:00		50ml ●	
CW-6601	light body fast		00:50	01:20		02:10		50ml ●	
CW-6511	regular body		01:00	02:00		03:00		50ml ●	
CW-6611	regular body fast		00:50	01:20		02:10		50ml ●	
CW-6775	PRECIOUS light body		01:00	02:20		03:20		50ml ●	
CW-6776	PRECIOUS regular body		01:00	02:20		03:20		50ml ●	
CW-6547	MonoBody		01:00	02:30		03:30		75ml ●	
CW-6482			01:00	02:30		03:30		380ml Refill ●	
CW-6520	heavy body		01:00	02:00		03:00		75ml ●	
CW-6487			01:00	02:00		03:00		380ml Refill ●	
CW-6620	heavy body fast		00:50	01:20		02:10		75ml ●	
CW-60019554	Black Edition		01:00	02:20		03:20		75ml ●	
CW-60019776	(heavy body)		01:00	02:20		03:20		380ml Refill ●	
CW-6535	putty super soft		01:05	02:00		03:05	600ml		
CW-6530	putty soft		01:10	02:00		03:10	600ml		
CW-6531	putty soft fast		01:05	01:50		02:55	600ml		

Accessories

Art. Nr.	Product		pcs
CW-60019663	Universal Mixing Tip (12x ●)		40
CW-60019694			100
CW-6552	Mixing Tip Turquoise (8x ●)		40
CW-6555	Oral Tips Yellow		100
CW-6162			Dynamic Mixing Tips Yellow (1x ●)
CW-6163	Fixation Ring		1

© COLTENE – www.coltene.com

Coltene/Whaledent Ltd.
The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF / UK
T +44 1444 235486
F +44 1444 870640
info.uk@coltene.com

Coltene/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten/Switzerland
T +41 71 757 5300
F +41 71 757 5301
info.ch@coltene.com

Coltene/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223/USA
T +1 330 916 8800
F +1 330 916 7077
info.us@coltene.com

 COLTENE