

For the ideal reproduction of function due to natural morphology.



## > FACTS & FEATURES

- What?**
- natural, detailed premium posterior teeth made of MRP composite
- Benefits?**
- **simple restoration of natural chewing function**  
thanks to occlusal surfaces with natural morphology
  - **intuitive centric finding and good occlusion**  
thanks to centric-supported occlusal design
  - **good shade reproduction with natural remaining dentition**  
thanks to excellent shade accuracy to the shade guide and broad shade selection
  - **reliable masking of tertiary structures**  
thanks to high (basal) tooth volume
  - **excellent durability**  
due to high load capacity and abrasion stability
- What for?**
- ideal solution for partial and hybrid dentures

# VITA PHYSIODENS®



## > BASIC INFORMATION

### Characteristics

#### Mould

- three-dimensional tooth geometry based on natural tooth dimension
- occlusal surfaces/cusp relief according to the natural, fully anatomical template

#### Special features

- high (basal) volume
- Maxillary cusp angle: approx.  $> 30^\circ$
- detailed fissure design

### Recommended indications

#### Ideal prosthetic solution for:

- partial dentures
- model casting, attachment and bar constructions
- Telescope and cone constructions

### Occlusal/setup concept

#### ★ Especially suited for:

BLP (BIO-Logical Prosthetic) according to Dr. E. End

### Available shades

#### • 27 x VITA SYSTEM 3D-MASTER® shades

0M3, 1M1–5M3

#### • 15 x VITA classical A1–D4® shades

A1–D4 (without B1)

### Available moulds/sets

- 6 x upper posteriors tooth moulds
- 6 x lower posterior tooth moulds

POSTERIOR

VITAPAN  
EXCELL

VITAPAN  
PLUS

VITAPAN  
LINGOFORM

VITA  
PHYSIODENS

VITA  
MIFT

CAD/CAM  
TEETH











CAD/CAM  
MATERIALS

CAD/CAM  
COMPONENTS



# VITA PHYSIODENS®

## > OVERVIEW OF SHAPES

Posterior						
<b>20E</b>	<b>21E*</b>	<b>22E*</b>	<b>23E*</b>	<b>24E*</b>	<b>25E</b>	
8.9	9.5	9.5	11.1	11.8	11.2	
						
						
8.3	8.7	9.0	9.8	9.8	9.5	

\*) Bleached shade 0M3 available exclusively for selected mould(s).

## > RECOMMENDED COMBINATIONS

### Combinations according to indications:

#### For partial and hybrid dentures

- VITA PHYSIODENS Posterior with VITA PHYSIODENS Anterior
- See combination table for VITA PHYSIODENS Anterior

## › ASSORTMENTS AND STORAGE



### VITA PHYSIODENS mould chart

Handy metal case with all moulds for tooth selection:

- VITA PHYSIODENS Anterior and Posterior



### VITA PHYSIODENS tooth storage

Posterior:

- X-small\*: 50 Sets / 400 Posteriors
- Small\*: 100 Sets / 800 Posteriors
- Medium\*: 200 Sets / 1,600 Posteriors
- Large\*: 300 Sets / 2,400 Posteriors

Available in VITA classical A1–D4® and VITA SYSTEM 3D-MASTER® shades.\*\*

\*) Note: You will receive the most frequently used moulds and shades with the packaged tooth storage. We are also happy to put together storage systems according to your individual specifications.

\*\*) see available shades

## › SYSTEM COMPONENTS

