


INSTI-DAM[®]



A QUICK WAY TO ISOLATE YOUR TREATMENT AREA & TAKE CONTROL OF AEROSOL CONTAMINANTS

Insti-Dam[®] is quick & easy to place, allowing clinicians to isolate their treatment area & take control of aerosol contaminants & aspiration hazards, thus creating a safer environment for both patients & clinicians.

INSTI-DAM®

Patent #7044738

Features

Pre-assembled with a built-in flexible frame

Allows radiographs to be taken by folding Insti-Dam® to the side

1/2" offset pre-punched hole

Makes positioning easy & helps eliminate tearing (additional holes can be punched)

Compact design

- 4-1/4" diameter; fits outside patients lips
- Easy to use
- Non-threatening and comfortable for patients


Medium gauge material

Disposable

No risk of cross contamination & no maintenance

Pro Tip

Pair with SE Cushion or Pink Petal® on a saliva ejector for extra patient comfort & overall saliva reduction!



Pre-punched hole

Pre-assembled on frame



Material can be adjusted on/around the frame once placed. (Relaxed Fit version ONLY)

Original Insti-Dam®

Ideal for anterior isolation

Insti-Dam® Original (Blue - 20pk)

- Made with polyisoprene
- Dimensions: 4-1/4" diameter

Insti-Dam® Latex (Natural - 20pk)

- Made with natural rubber latex
- Dimensions: 4-1/4" diameter

Relaxed-Fit Insti-Dam®

Ideal for posterior isolation

Insti-Dam® Relaxed Fit (Purple - 20pk)

- Made with silicone
- The extra "relaxed" material in the Relaxed-Fit Insti-Dam® makes for quick placement with no stretching, tearing, or pulling.
- Dimensions: 4-7/8" x 4-7/8" x 6-1/8"



Because Time Is Everything

www.zirc.com